

## Amine-PEG3-OH

<http://www.lumiprobe.com/p/amine-peg3-oh-6338-55-2>

Bifunctional PEG3 derivative with amine and hydroxyl groups. The molecule contains an amine group that is reactive towards many electrophilic reagents like activated esters and terminal hydroxyl. This is a hydrophilic compound miscible with water.

---

### General properties

Appearance: colorless liquid

Molecular weight: 149.19

CAS number: 6338-55-2

Molecular formula:  $C_6H_{15}NO_3$

IUPAC name: 2-[2-(2-aminoethoxy)ethoxy]ethanol

Solubility: miscible with water, polar organic solvents

Quality control: 95+% GCMS,  $^1H$  NMR

Storage conditions: Storage: 24 months after receipt at  $-20^\circ C$  in the dark. Transportation: at room temperature for up to 3 weeks.

Legal statement: This Product is offered and sold for research purposes only. It has not been tested for safety and efficacy in food, drug, medical device, cosmetic, commercial or any other use. Supply does not express or imply authorization to use for any other purpose, including, without limitation, in vitro diagnostic purposes, in the manufacture of food or pharmaceutical products, in medical devices or in cosmetic products.