

## **QuReader® 1 Green**

<http://www.lumiprobe.com/p/fluorometer-quireader-1-well-green>

The QuReader® 1 Green fluorometer contains one fluorescence channel and enables dsDNA, ssDNA, RNA, and protein quantitation using fluorescence measuring in green channel. Total compatibility with the QuDye® assays ensure highly sensitive determination of nucleic acids and protein concentrations in the lower pg/uL range.

---

### **General properties**

**Legal statement:** This Product is offered and sold for research purposes only. It has not been tested for safety and efficacy in food, drug, medical device, cosmetic, commercial or any other use. Supply does not express or imply authorization to use for any other purpose, including, without limitation, in vitro diagnostic purposes, in the manufacture of food or pharmaceutical products, in medical devices or in cosmetic products.