

## **QuReader® 8 Green**

<http://www.lumiprobe.com/p/fluorometer-quireader-8-well-green>

The QuReader® 8 Green fluorometer is a desk fluorometer with a touch screen designed for dsDNA, ssDNA, RNA, and protein quantitation. The QuReader® fluorometer enables the fluorescence measuring for up to eight samples using green channel. Total compatibility with the QuDye® assays ensure highly sensitive determination of nucleic acid and protein concentrations in the lower pg/uL range.

---

### **General properties**

**Legal statement:** This Product is offered and sold for research purposes only. It has not been tested for safety and efficacy in food, drug, medical device, cosmetic, commercial or any other use. Supply does not express or imply authorization to use for any other purpose, including, without limitation, in vitro diagnostic purposes, in the manufacture of food or pharmaceutical products, in medical devices or in cosmetic products.