


(C5) Isovaleryl-L-carnitine-(trimethyl-d9)

<http://www.lumiprobe.com/p/isovaleryl-l-carnitine-trimethyl-d9>

Isovaleryl-L-carnitine-(N,N,N-trimethyl-d9) is a standard for the quantitative determination by MS/MS of C5 isovalerylcarnitine, C5:1 tiglicarnitine, C5OH/C4DC ratio 3-hydroxyisovalerylcarnitine/methylmalonocarnitine.



Structure of (C5) Isovaleryl-L-carnitine-(trimethyl-d9)

General properties

- Appearance: white solid
- Molecular weight: 290.83
- CAS number: 203806-03-5
- Molecular formula: $C_{12}H_{15}D_9ClNO_4$
- Solubility: DMF, DMSO, ethanol
- Quality control: NMR 1H and HPLC-MS (95+ %, D: 98+ %)
- Storage conditions: 24 months after receipt at -20°C in the dark. Transportation: at room temperature for up to 3 weeks. Desiccate.
- Legal statement: This Product is offered and sold for research purposes only. It has not been tested for safety and efficacy in food, drug, medical device, cosmetic, commercial or any other use. Supply does not express or imply authorization to use for any other purpose, including, without limitation, in vitro diagnostic purposes, in the manufacture of food or pharmaceutical products, in medical devices or in cosmetic products.