


## L-Ornithine-d2

<http://www.lumiprobe.com/p/l-ornithine-d2>

L-Ornithine-d2 is the deuterium labeled amino acid L-Ornithine. Stable heavy isotopes such as hydrogen are used for labeling for quantitation in pharmacokinetic and metabolic application as tracers.



**Structure of L-Ornithine-d2**

---

### General properties

Appearance: white powder

Molecular weight: 134.18

Molecular formula: C<sub>5</sub>H<sub>10</sub>D<sub>2</sub>N<sub>2</sub>O<sub>2</sub>

Solubility: in water

Quality control: NMR <sup>1</sup>H and HPLC-MS (95+ %, D: 98+ %)

Storage conditions: 24 months after receipt at -20°C in the dark. Transportation: at room temperature for up to 3 weeks. Desiccate.

Legal statement: This Product is offered and sold for research purposes only. It has not been tested for safety and efficacy in food, drug, medical device, cosmetic, commercial or any other use. Supply does not express or imply authorization to use for any other purpose, including, without limitation, in vitro diagnostic purposes, in the manufacture of food or pharmaceutical products, in medical devices or in cosmetic products.