


L-Phenylalanine (ring-d5)

<http://www.lumiprobe.com/p/l-phenylalanine-ring-d5>

L-Phenylalanine-d5 is a stable isotope of phenylalanine can be used to study the metabolic effects and proteins expression by the analytical methods using NMR.



Structure L-Phenylalanine (ring-d5)

General properties

Appearance: white powder

Molecular weight: 170.22

CAS number: 56253-90-8

Molecular formula: C₉H₆D₅NO₂

Solubility: in water

Quality control: NMR ¹H and HPLC-MS (95+ %, D: 98+ %)

Storage conditions: 24 months after receipt at -20°C in the dark. Transportation: at room temperature for up to 3 weeks. Desiccate.

Legal statement: This Product is offered and sold for research purposes only. It has not been tested for safety and efficacy in food, drug, medical device, cosmetic, commercial or any other use. Supply does not express or imply authorization to use for any other purpose, including, without limitation, in vitro diagnostic purposes, in the manufacture of food or pharmaceutical products, in medical devices or in cosmetic products.